

AL-2000 Data format

SOH		I D										DATE										T O M E Y			SPACE			TYPE					PATIENT NAME										AGE		SEX	PHYSICIAN NAME										C R	
0		0 0 0 0 0 1 2 3 4 5 6										2 0 0 1 0 8 3 1 1 1 3 5 2 0										T O M E Y			_ _ _ _ _			A L - 2 0 0 0					M R . A A A										3 5		M	D R . T O M E Y										C R	
STX	B M	C	R	R: right	L: left	E Y E T Y P E										AverageVelocity			LensVelocity			C		R																																	
STX	H	R	R	ACDVelocity	BioVerocity																																																				
STX	V	R	R	Axial Length	ACD	LENS																																																			
STX	L	R	R	K1	K2																																																				
STX	K	R	R	Desired Ref																																																					
STX	D	R	R	IOL Thickness																																																					
STX	T	R	R	IOL formula																																																					
STX	F	R	R	IOL Company	IOL model																																																				
STX	I	R	R	LensConstant1	LensConstant2	LensConstant3	Power										REF																																								
STX	R	R	R	IOL Company	IOL model																																																				
STX	I	R	R	LensConstant1	LensConstant2	LensConstant3	Power										REF																																								
STX	R	R	R	IOL Company	IOL model																																																				
STX	I	R	R	LensConstant1	LensConstant2	LensConstant3	Power										REF																																								
STX	R	R	R	E Y E T Y P E	AverageVelocity										LensVelocity			C		R																																					
STX	H	L	R	ACDVelocity	BioVerocity																																																				
STX	V	L	R	Axial Length	ACD	LENS																																																			
STX	L	L	R	K1	K2																																																				

